

DAB Substrate Kit

Catalogue No.: IHC0005

(1 Kit(20X) for 1200 assays)

FineTest DAB Substrate (3, 3'-diaminobenzidine) produces a brown reaction product in the presence of peroxidase (HRP) enzyme. Substrate Kit contains all of the reagents necessary to prepare. These reagents are supplied in convenient dropper bottles.

REAGENTS:

Reagents A	DAB 20×concentrated solution	3ml	-20°C
Reagents B	H ₂ O ₂ 20×concentrated solution	3ml	4-8°C
Reagents C	TBS 20× concentrated solution	3ml	4-8°C

STORAGE:

- Store reagents in original bottles at the required conditions (one year)
- Avoid storing reagents or working solution in strong direct light

PREPARATION OF SUBSTRATE WORKING SOLUTION:

- To 1.0 ml of distilled water
- Add 2 drops (approximately 50 µl) of Reagents C
- Add 2 drops (approximately 50 µl) of Reagents A
- Add 2 drops (approximately 50 µl) of Reagents B
- Mix well before use

INSTRUCTIONS FOR USE:

For Tissues or Cells

After incubation with a peroxidase (HRP) detection system, rinse well. Incubate with the substrate working solution at room temperature for 2-10 minutes. Optimal development times should be determined by the investigator.

Wash for 5 minutes in water.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com

NOTES:

- We recommend using glass-distilled water in the preparation of substrate buffer. Deionized water may contain inhibitors of the peroxidase reaction. Solutions containing sodium azide or other inhibitors of peroxidase activity should not be used in diluting the peroxidase substrate.
- DAB are carcinogens. Appropriate care should be exercised when using these reagents including gloves, eye protection, lab coats, and good laboratory procedures.

Wuhan Fine Biotech Co., Ltd.

B9 Bld, High-Tech Medical Devices Park, No. 818 Gaoxin Ave. East Lake High-Tech Development Zone. Wuhan, Hubei, China(430206)

Tel : (0086)027-87384275

Fax: (0086)027-87800889

www.fn-test.com